


# REBOA Instructions

This page is intended to serve as a quick reference for easy access to the REBOA kit instructions. The information and images are directly from the insertion instructions, and were obtained from <http://prytime.com/wp-content/uploads/2017/08/ER-REBOA-Catheter-Quick-Reference-Guide-wall-poster.pdf>


## The ER-REBOA™ Catheter Quick Reference Guide 6 REBOA Steps: ME-FIIS (Pronounced ME-FIZZ)

### Get Access Early



Obtain access using standard techniques

### 1. Measure



Placement depth<sup>1,2,3,4,5,6</sup>


- Zone 1: Approximately 46 cm
- Zone 3: Approximately 28 cm

### 2. Empty



Flush & deflate balloon

- Ensure balloon is fully deflated
- Hold vacuum for **5 seconds**
- Close stopcock with vacuum held



Advance & twist peel-away to cover P-tip®

- Ensure the balloon and P-tip® are captured


### 3. Flush



Attach & flush arterial line


- Use standard techniques
- Ensure all air is purged

### 4. Insert



Insert sheath into valve

- Approximately 5 mm
- Insert into the common femoral artery



Advance catheter into vessel

- Hold orange sheath
- Advance blue Catheter
- Remove sheath after balloon passes valve



Position catheter

If available, use conventional x-ray or fluoroscopy to confirm position using radiopaque markers


## 5. Inflate<sup>1,2,3,4,5,6</sup>

Inflation Volume

**Zone 1** Start with 8 cc

**Zone 3** Start with 2 cc


### Start small then check



#### Monitor arterial waveform feedback

- Look for change in blood pressure above balloon
- Use other standard techniques

## 6. Secure



#### Secure Catheter close to the introducer sheath

## Provide Definitive Treatment



#### Provide definitive hemorrhage control

- Mark time of inflation
- The clock is ticking!
- Move quickly to definitive control

## Remove



#### Fully deflate balloon

- Hold vacuum for **5 seconds**
- Close stopcock with vacuum held



#### Remove catheter

- Corkscrew twist the catheter to facilitate removal
- If necessary, remove catheter and introducer sheath as a unit

## Caution



#### Check for full and equal pulse in each leg using your standard technique



This instruction is not a replacement for the instruction for use (IFU). The ER-REBOA™ Catheter IFU should be read in its entirety before using the device.

1. Joint Trauma System Clinical Practice Guideline (JTS CPQ) REBOA for Hemorrhagic Shock (CPQ-2). 2018.
2. Papp A, Barak A, Pore N, Conon F, Lundberg PM, Cellat J, David J, Volpic G. Fixed-Diameter Voids for Balloon Placement During Fluoroscopically-Guided Endovascular Balloon Occlusion of the Aorta in a Civilian Population. *JVIR*. 2018;20(1):14-16.
3. Lindebar M, Traub K, Helmke T, Rasmussen TE, Smith J, Wendeborg R, Grato D, Demetriadou G. Emergent non-image-guided revascularization balloon occlusion of the aorta (REBOA) for major trauma: A cadaver-based study. *J Trauma Acute Care Surg*. 2018;84(2):242-7.
4. MacTaggart JM, Poulsen WG, Alomar JJ, Saxe A, Thorson K, Phillips M, Desjardis AG. Randomized trial: fluoroscopic roadmaps to improve accurate device delivery for fluoroscopically-guided endovascular balloon occlusion of the aorta. *J Trauma Acute Care Surg*. 2018;84(5):941-6.
5. Morrison JJ, Scannard A, Midwinter MJ, Sharon CJ, Elison JL, Rasmussen TE. Preoperative evaluation of the correlation between torso height and aortic anatomy in reboas of a fluoroscopically-guided balloon occlusion system. *Surgery*. 2014 Jun;156(2):244-9.
6. Scannard A, Morrison J, Sharon CJ, Elison JL, Rasmussen TE. Ultrasonographic analysis of torso aortic anatomy with implications for revascularization balloon occlusion. *J Trauma Acute Care Surg*. 2015 Aug;78(2 Suppl):2516-72.

ADV-006 | Revision E

Revision #3

Created 28 August 2023 18:12:42 by Abby Josef

Updated 12 February 2026 15:49:39 by Abby Josef